



# Stabilis



## Metronidazole benzoate



### Stabilité des préparations

1400 mg Metronidazole benzoate		Propylène Glycol 33,3 mg SyrSpend SF® >> 20 ml	2-8°C		730			2250
1400 mg Metronidazole benzoate		Propylène Glycol 33,3 g SyrSpend SF® >> 20 ml	25°C		365			2250





### Bibliographie

Type	Source
2250 Revue	Vu N.T, Aloumanis V, Ben MJ, Kupiec T.C, Patterson E K, Radke J, Erickson M A III, Schneider G. Stability of Metronidazole Benzoate in SyrSpend SF One-Step Suspension System. Int J Pharm Compound 2008 ; 12, 6: 558-564.



# Dictionnaire

 Antibiotique	 Solution buvable
 Stabilité des préparations	 Contenant
 Origine	 Excipient
 Température	 Conservation
 Durée de stabilité	 Biosimilaire
 Données conflictuelles	 Bibliographie
 Flacon plastique	 Poudre
 A l'abri de la lumière	 Jour
 Bibliographie	 Dictionnaire