














































# Stabilis



								
Aztreonam	0,4 mg/ml	PVC	Nutri neal®	25°C		3		4456
Aztreonam	0,4 mg/ml	PVC	Nutri neal®	37°C		1		4456
Aztreonam	0,4 mg/ml	PVC	Nutri neal®	6°C		14		4456
Ciprofloxacin lactate	47,6 mg/ml	PVC	Nutri neal®	25°C		7		4460
Ciprofloxacin lactate	47,6 mg/ml	PVC	Nutri neal®	37°C		24		4460
Ciprofloxacin lactate	47,6 mg/ml	PVC	Nutri neal®	6°C		14		4460
Fosfomicin	1,58 mg/ml	PVC	Nutri neal®	25°C		14		4458
Fosfomicin	1,58 mg/ml	PVC	Nutri neal®	37°C		24		4458
Fosfomicin	1,58 mg/ml	PVC	Nutri neal®	6°C		14		4458
Linezolid	0,21 mg/ml	PVC	Nutri neal®	37°C		24		4176



# Dictionnaire

 Listes récapitulatives	 Stabilité des solutions
 Solvant	 Nutrineal®
 Molécule	 Concentration
 Contenant	 Température
 Conservation	 Durée de stabilité
 Bibliographie	 Polyvinyl chlorure
 A l'abri de la lumière	 Jour
 Heure	 Non précisée
 Dictionnaire	