







































































Stabilis



|  |  |  |  |  |  |  |  | |
|---|---|---|---|--|---|---|---|-----|
| Cytarabine | 3 mg/ml |  | Ei B | 23°C |  | 48 |  | 31 |
| Cytarabine | 3 mg/ml |  | Ei B | 4°C |  | 48 |  | 31 |
| Cytarabine | 5 mg/ml |  | Ei B | 25°C |  | 7 |  | 739 |
| Cytarabine | 3 mg/ml |  | Ei B | 23°C |  | 48 |  | 31 |
| Cytarabine | 3 mg/ml |  | Ei B | 4°C |  | 48 |  | 31 |
| Hydrocortisone sodium succinate | 1 mg/ml |  | Ei B | 25°C |  | 24 |  | 739 |
| Hydrocortisone sodium succinate | 2 mg/ml |  | Ei B | 23°C |  | 24 |  | 31 |
| Hydrocortisone sodium succinate | 2 mg/ml |  | Ei B | 4°C |  | 48 |  | 31 |
| Hydrocortisone sodium succinate | 2 mg/ml |  | Ei B | 23°C |  | 24 |  | 31 |
| Hydrocortisone sodium succinate | 2 mg/ml |  | Ei B | 4°C |  | 48 |  | 31 |
| Methotrexate sodium | 2 mg/ml |  | Ei B | 23°C |  | 48 |  | 31 |
| Methotrexate sodium | 2 mg/ml |  | Ei B | 4°C |  | 48 |  | 31 |
| Methotrexate sodium | 2,5 mg/ml |  | Ei B | 25°C |  | 7 |  | 739 |
| Methotrexate sodium | 2 mg/ml |  | Ei B | 23°C |  | 48 |  | 31 |
| Methotrexate sodium | 2 mg/ml |  | Ei B | 4°C |  | 48 |  | 31 |



Dictionnaire

| | |
|--|---|
|  Listes récapitulatives |  Stabilité des solutions |
|  Solvant | El B Solution d'Elliott B |
|  Molécule |  Concentration |
|  Contenant |  Température |
|  Conservation |  Durée de stabilité |
|  Bibliographie |  Verre |
| ? Non précisée |  Heure |
|  A l'abri de la lumière |  Lumière |
|  Jour |  Seringue polypropylène |
|  Dictionnaire | |