



















































Stabilis



|  |  |  |  |  |  |  |  | |
|---|---|---|---|--|---|---|---|------|
| Linezolid | 0,15 -> 0,6 mg/ml |  |  | 25°C |  | 7 |  | 2163 |
| Linezolid | 0,15 -> 0,6 mg/ml |  |  | 37°C |  | 7 |  | 2163 |
| Linezolid | 0,15 -> 0,6 mg/ml |  |  | 4°C |  | 7 |  | 2163 |
| Piperacillin sodium / tazobactam | 0.2 / 0.025 mg/ml |  |  | 23°C |  | 7 |  | 2451 |
| Piperacillin sodium / tazobactam | 0.2 / 0.025 mg/ml |  |  | 37°C |  | 24 |  | 2451 |
| Piperacillin sodium / tazobactam | 0.2 / 0.025 mg/ml |  |  | 4°C |  | 14 |  | 2451 |



Dictionnaire

| | |
|--|---|
|  Listes récapitulatives |  Stabilité des solutions |
|  Solvant |  Dianéal® PD2 4.25% glucose (Baxter) |
|  Molécule |  Concentration |
|  Contenant |  Température |
|  Conservation |  Durée de stabilité |
|  Bibliographie |  Polyvinyl chlorure |
|  A l'abri de la lumière |  Jour |
|  Non précisé |  Non précisée |
|  Heure |  Dictionnaire |