

















































# Stabilis



								
Amoxicillin sodium	0.125 mg/ml	PVC		25°C	?	72		3777
Amoxicillin sodium	0.125 mg/ml	PVC		37°C	?	24		3777
Amoxicillin sodium	0.125 mg/ml	PVC		4°C	?	14		3777
Ampicillin sodium	0.125 mg/ml	PVC		25°C	?	72		3777
Ampicillin sodium	0.125 mg/ml	PVC		37°C	?	24		3777
Ampicillin sodium	0.125 mg/ml	PVC		4°C	?	14		3777
Insulin	0,004 UI/ml	PVC		25°C	?	24		2164
Insulin	0,01 UI/ml	PVC		25°C	?	24		2164
Insulin	0,02 UI/ml	PVC		25°C	?	24		2164
Insulin	0,04 UI/ml	PVC		25°C	?	24		2164
Oxaliplatin	0,1 mg/ml			42°C	?	0.5		2791



# Dictionnaire

 Listes récapitulatives	 Stabilité des solutions
 Solvant	 Dianéal® PD4 1.36% glucose (Baxter)
 Molécule	 Concentration
 Contenant	 Température
 Conservation	 Durée de stabilité
 Bibliographie	 Polyvinyl chlorure
 Non précisée	 Heure
 Jour	 Verre
 Dictionnaire	